

\*\*\*When Starting @100F\*\*\*

# MilkWorks

by Dairy Tech, Inc.™



## Thaw and Warm Times



PERFECT UDDER BAGS # of bags	time in MINUTES	*2qt CALF BOTTLES # of bottles
1 to 2	25	
2 to 3	27	
3 to 4	30	
4 to 5	35	
5+	40	
	45	
	50	1 to 3
	55	4 to 6
	60	7 to 10

Always double check the temperature of the product before feeding calves. Cold colostrum has poor absorption. Hot colostrum may burn the esophagus severely.

*For more information visit [www.Milk-Works.info](http://www.Milk-Works.info)*

\*3qt bottles will require more time